


Fast Track Media Fine-tuning Service

| Key features

- Tailored media and/or supplement customization for enhancing cell culture performance and process efficiency (non-GMP media prototypes).
- Flexible production sizes according to customer requirements.
- Exceptionally swift prototype production and delivery, typically completed within 10 business days.

| Overall Workflow



| We provide:

- Powder media/supplement customized and optimized to your cell line (including any specific requests for formulation adjustment)
- Media prototype delivered in lot size of up to 10kg, in aluminum pouches or pails, according to weight

| Quality control testing

- Result of Analysis (RoA): Provided to ensure the quality and performance of the custom media.
- Testing Parameters: Include appearance, solubility, pH, and osmolality.
- Optional Testing: TAMC/TYMC, endotoxin, mycoplasma, and cell culture testing are available but may require additional lead-time.

| Requirements from customer

- Spent Media Analysis (SMA) Data: Customers are encouraged to provide detailed results from their initial trials with our catalog media, including viable cell density (VCD), viability, and metabolite concentrations during culture. Protein titer data is optional.
- Spent Media Samples: Alternatively, customers can send spent media samples from their trials to CELLiST for comprehensive analysis. This allows us to perform SMA, measuring all metabolites and amino acid concentration changes during the culture.
- Additional Information: Any specific requirements or constraints related to the customer's cell line and process to tailor the media prototype effectively.

Sales & Technical Support

For quoting, ordering, product sample request

CELLiST™ Sales Team

Tel: +82-032-210-2600

Mail: sales@ajinomotocellist.com

Global CELLiST Solution Centers



CELLiST | Solution Center Korea

70, Songdogwahak-ro, Yeonsu-gu, Incheon, Republic of Korea

Tel: +82-32-210-2695 Fax: +82-32-210-2604

Mail: yaron.silberberg.dk2@asv.ajinomoto.com

CELLiST | Solution Center India

2nd Floor, Block 05, TIMES CORPORATE PARK Thaltej, Ahmedabad, Gujarat 380058, India

Tel: +91-87359-10813

Mail: parshadg.hirapara.ut3@asv.ajinomoto.com

CELLiST | Solution Center China

Room 501, Building 9, NO.185 Jiangtian East Road, Songjiang District, Shanghai, China

Tel: +86-21-6219-2705

Mail: mingqi_fan@ajinomoto.com.cn

Eat Well, Live Well.

